



This table outlines known interferences for the assays on the Siemens epoc BGEM test card. As with any laboratory test, other unknown interferences may affect patient results. Consult a Clinical Biochemist if assay interference is suspected but not listed below.

Analyte	Interfering Substances	Effect on Results
pH	Benzalkonium salts used in coatings for in-dwelling lines. Lines must be properly flushed prior to sampling.	Falsely decreases pH
	Exposure to air	Falsely increases pH
pCO ₂	Bromide salts	Increases pCO ₂ by 0.2 mmHg/mM Br
	Exposure to air	Falsely decreases pCO ₂
pO ₂	Metronidazole	Increases pO ₂ by 4 mmHg at 100 µM Metronidazole (therapeutic range)
	Exposure to air	Falsely increases pO ₂
Sodium	High total protein (e.g. multiple myeloma)	Decreases Na ⁺ by 1.3 mM for every additional 10 g/L protein above 80 g/L
	β-hydroxybutyrate	20 mmol/L βHB will decrease Na ⁺ by 3 mmol/L. Results above 18 mmol/L have not been reported in Manitoba.
	Bromide	16 mmol/L Bromide will increase Na ⁺ by 5 mmol/L
	Benzalkonium salts used in coatings for in-dwelling lines. Lines must be properly flushed prior to sampling.	Causes significant elevations in Na ⁺
	Teriflunomide	Causes decrease in Na ⁺ at supra-therapeutic levels (>720 µmol/L)
Potassium	Benzalkonium salts used in coatings for in-dwelling lines. Lines must be properly flushed prior to sampling.	Causes significant increase in K ⁺
	Excessive pulling on syringe plunger; use of small-bore needles; aggressive sample mixing—all cause hemolysis <i>in vitro</i>	Falsely increases K ⁺
Ionized Calcium	β-hydroxybutyrate (<i>extremely elevated</i>)	20 mmol/L βHB will decrease iCa by 0.038 mmol/L. Results above 18 mmol/L have not been reported in Manitoba.
	Salicylates	4.3 mmol/L salicylate (supratherapeutic levels) will decrease iCa by 0.06 mmol/L
	Bromide	10 mmol/L Bromide will increase iCa by 0.05 mmol/L
	Sodium perchlorate	1 mmol/L perchlorate will decrease iCa by 0.23 mmol/L
	Leflunomide	Causes decrease in iCa at therapeutic levels
	Teriflunomide	Causes decrease in iCa at therapeutic levels



Analyte	Interfering Substances	Effect on Results
Chloride	β-hydroxybutyrate	βHB above 8.3 mmol/L will increase chloride by 0.37 mmol/L for every additional 1 mmol/L βHB. βHB results this high are only observed in severe diabetic ketoacidosis.
	Bromide	Bromide will increase Cl ⁻ by up to 2.75 mmol/L per mmol/L Bromide
	Citrate	Citrate may increase Cl ⁻ by up to 1.13 mmol/L per mmol/L citrate
	N-acetylcysteine (NAC)	Supratherapeutic concentrations may decrease Cl ⁻ by up to 5.55 mmol/L per mmol/L NAC
	Salicylates	Supratherapeutic concentrations may decrease Cl ⁻ by 2.94 mmol/L for every 1 mmol/L salicylate above 1.67.
	Thiocyanate (SCN ⁻)	Can increase Cl ⁻ by up to 5.62 mmol/L for every mmol/L Thiocyanate
Total CO₂	Exposure to air	Falsely decreases TCO ₂
	Benzalkonium salts used in coatings for in-dwelling lines. Lines must be properly flushed prior to sampling.	Causes significant decrease in TCO ₂
	Citrate	Can increase TCO ₂ by up to 0.24 mmol/L per mmol/L citrate
	EDTA	Can increase TCO ₂ by up to 0.57 mmol/L per mmol/L EDTA
	N-acetylcysteine (NAC)	Supratherapeutic concentrations may increase TCO ₂ by up to 0.54 mmol/L per mmol/L NAC
Glucose	Capillary collection	Typically, 0.1-0.3 mmol/L higher than venous or arterial glucose
	Potassium Iodide	May falsely decrease glucose concentration by up to 0.01 mmol/L per μmol/L KI
	Bromide	May falsely decrease glucose by 0.013 mmol/L per mmol/L Bromide
	N-acetylcysteine (NAC)	Supratherapeutic concentrations may decrease glucose by up to 0.4 mmol/L per mmol/L NAC.
	Thiocyanate (SCN ⁻)	May decrease glucose by up to 1.7%
	Uric Acid	Uric acid concentrations above 700 μmol/L may decrease glucose concentrations by 0.2 mmol/L for every μmol/L uric acid over this threshold

Analyte	Interfering Substances	Effect on Results
Glucose (cont'd)	Mannose, Xylose	Both may falsely increase glucose concentration
Lactate	Acetaminophen	May falsely increase lactate in overdose situations. Expected increase of 0.3 mmol/L for every 150 mg/L acetaminophen above 120.
	Iodide	May falsely decrease lactate readings by up to 1.0 mmol/L
	Bromide	May falsely decrease lactate concentration by up to 0.15 mmol/L for every mmol/L Br.
	Thiocyanate (SCN ⁻)	May falsely decrease lactate concentration by up to 0.97 mmol/L for every mmol/L KSCN
	N-acetylcysteine (NAC)	Supratherapeutic concentrations may decrease lactate by up to 0.96 mmol/L per mmol/L NAC.
Urea	Benzalkonium salts used in coatings for in-dwelling lines. Lines must be properly flushed prior to sampling.	Causes significant increase in urea
	Citrate	Can decrease urea readings by up to 0.1 mmol/L per mmol/L citrate
	EDTA	Can decrease urea readings by up to 0.15 mmol/L per mmol/L EDTA
	Glutathione	Long-term oral glutathione or IV glutathione administration may falsely increase urea concentrations by 0.68 mmol/L per mmol/L GSH. Short-term oral intake should not affect results.
	β -hydroxybutyrate (<i>extremely elevated</i>)	β HB above 17 mmol/L will increase urea by 0.05 mmol/L for every additional 1 mmol/L β HB. Results above 18 mmol/L have not been reported in Manitoba.
	Intravenous hydroxyurea	Plasma concentrations of hydroxyurea above 1.3 mmol/L may increase Urea by 0.6 mmol/L per mmol/L hydroxyurea above this threshold.
	N-acetylcysteine (NAC)	Supratherapeutic concentrations (>9 mmol/L) may decrease urea by up to 0.04 mmol/L per mmol/L NAC.
	Sodium Thiosulfate (Na ₂ S ₂ O ₃)	Concentrations above 4.1 mmol/L can decrease urea concentration by up to 0.15 mmol/L per mmol/L sodium thiosulfate.



Analyte	Interfering Substances	Effect on Results
Creatinine	Creatine	Will falsely increase creatinine concentration if consumed within 48 hours of testing
	Iodide	Plasma concentrations above 0.45 mmol/L may decrease creatinine concentrations by up to 44 µmol/L for every mmol/L iodide.
	N-acetylcysteine (NAC)	Therapeutic concentrations (>0.4 mmol/L) may decrease creatinine concentration by up to 64 µmol/L per mmol/L NAC.
	Glycolic acid (metabolite of ethylene glycol)	Plasma concentrations above 1.69 mmol/L can decrease creatinine by up to 5%.
Hematocrit	Abnormally high lipids	Falsely increase hematocrit results
	Abnormally high white blood cell count	May falsely increase hematocrit results
	Total protein content	Normal plasma total protein content is 46 – 80 g/L accounting for all age groups. A 10 g/L increase in total protein can falsely increase hematocrit by 0.01 L/L. Similarly, a 10 g/L decrease in total protein can falsely decrease hematocrit by approximately 0.01 L/L.