

For PDF Fillable Requisitions, the following applies:

1. The form shall be completed using a Digital Health assigned computer.
2. Absolutely no personal health information shall be electronically saved on a computer.
3. The completed form shall not be shared electronically. If you reasonably believe that e-mailing the information is the only available method of communication or the only way to send the information then you must adhere to the Privacy guideline titled "E-mailing Personal Health Information".
4. All forms must be completed in their entirety, e.g. if a staff member has only completed half of a form they cannot save their work and then come back to complete it at a later date.
5. Once the personal health information has been recorded onto the form, it is to be printed immediately, deleted (not saved) from the computer, and then stored securely inside the client (paper) health record or scanned into an electronic record.
6. Do not print unnecessary duplicate copies of the form.
7. Regular audits of the Digital Health assigned computer shall be undertaken to ensure that no personal health information is being duplicated and saved.



DIAGNOSTIC SERVICES  
OF MANITOBA

SERVICES DE DIAGNOSTIC  
DU MANITOBA

ADDRESSOGRAPH PREFERRED  
or PRINT ONLY:

LOCATION:  
WARD

PATIENT NAME:  
LAST, FIRST

DATE OF BIRTH:  
DD/MM/YYYY

**BIOCHEMISTRY TEST REQUISITION**

SEX  F  M

MB PHIN:  
(Specify province if different)

PHYSICIAN: (PRINT)  
LAST, FIRST

PHYSICIAN BILLING CODE:

ORDERING PROFESSIONAL:  
(If different from physician)

Collection Date: \_\_\_\_\_ Time: \_\_\_\_\_  
d/m/y

**RASBURICASE PROTOCOL – URIC ACID**

**Specimen Collection**

**Adult: 1 – 4.5 mL Lithium Heparin PST Tube (\*PLACE PROMPTLY ON ICE)**

**Pediatric: 1 – 0.6 mL Lithium Heparin Microtainer PST Tube (\*PLACE PROMPTLY ON ICE)**

***KEEP ON ICE & TRANSPORT IMMEDIATELY TO BIOCHEMISTRY LABORATORY***

**Test Name**

**Test Code**

Rasburicase Protocol - Uric Acid

RPUA

**PROCESSING**

- Sample **MUST** be centrifuged at 4°C, keep ON ICE until analysis.
- Perform analysis without delay
- Program manually at the instrument for uric acid
- Enter the result on the RPUA format

**\*\*\*Caution: Rasburicase protocol – uric acid is for patients being treated with the enzyme Rasburicase. This enzyme breaks down uric acid ex-vivo. Thus, samples for rasburicase protocol must be kept at 2-8 °C and analyzed promptly.\*\*\***