

	Approx. % cross reactivity
3,4-Methylenedioxymethamphetamine (MDMA)	197
para-Methoxyamethamphetamine (PMMA)	145
para-Methoxyamphetamine (PMA)	110
3,4-Methylenedioxyamphetamine (MDA)	130
d-Amphetamine	102
d-Methamphetamine	100
N-methyl-1-(1,3-benzodioxol-5-yl)-2-butanamine (MDBD)	85
3,4-Methylenedioxy-N-ethylamphetamine (MDEA)	64
3,4-Methylenedioxyphenyl-2-butanamine (BDB)	41
Trazodone metabolite: mCPP	18
1-Methyl-3-phenylpropylamine	20
l-Methamphetamine	11
l-Amphetamine	4
Dimethylamylamine	0.99