

	Approx. % cross reactivity
Hydroxymethadone	9.1
Vortioxetine	4.1
Lu AA34443	11
Cyamemazine	3.5
Methotrimeprazine (Levomepromazine)	3.4
Chlorpromazine	1.2
Thiothixene	0.8
Clomipramine	0.2
Promazine	0.2
Thioridazine	0.2
Chlorprothixene	0.2
l- $\alpha$ -methadol	0.1
Promethazine	0.1
l- $\alpha$ -acetylmethadol (LAAM)	0.1
Trimipramine	0.1